



Twenty-Nine Palms Band of Mission Indians

29 Palms Laboratory

Laboratory Analysis Price List

2005-01-07

Pesticide Testing

Organochlorine (OC)	-	EPA Method 8081A	-	\$ 100.00
Organophosphate (OP)	-	EPA Method 8141A	-	\$ 150.00
Methyl Carbamate (MC)	-	EPA Method 8318	-	\$ 200.00
Multi-Screen (OC,OP, MC)	-	8081A/8141A/8318	-	\$ 400.00

Microbiology

Heterotrophic Plate Count	-	SM 9215 B	-	\$ 15.00
Total and Fecal Coliforms (P/A)	-	SM 9223 B	-	\$ 15.00
Total and Fecal Coliforms (MPN)	-	SM 9221	-	\$ 40.00
Total and Fecal Coliforms (MPN)	-	SM 9223 B	-	\$ 28.00

Metals

Metals Panel	-	EPA Method 6010	-	\$ 125.00
Aluminum (Al)	Lead (Pb)			
Antimony (Sb)	Magnesium (Mg)			
Arsenic (As)	Manganese (Mn)			
Barium (Ba)	Nickel (Ni)			
Beryllium (Be)	Potassium (K)			
Cadmium (Ca)	Selenium (Se)			
Calcium (Cd)	Silver (Ag)			
Chromium (Cr)	Sodium (Na)			
Cobalt (Co)	Thallium (Tl)			
Copper (Cu)	Vanadium (V)			
Iron (Fe)	Zinc (Zn)			

Plus an additional \$15.00 digestion fee per sample

Individual Metals - \$ 10.00 per metal with a per sample digestion fee



Twenty-Nine Palms Band of Mission Indians

29 Palms Laboratory

Nutrients and Minerals

Chloride	-	EPA Method 325.3	-	\$ 10.00
Fluoride	-	EPA Method 340.1	-	\$ 15.00
Nitrate	-	EPA Method 352.1	-	\$ 20.00
Phosphate, ortho	-	EPA Method 365.3	-	\$ 20.00
Phosphorus, total	-	EPA Method 365.3	-	\$ 25.00
Sulfate	-	EPA Method 375.4	-	\$ 20.00

Anions

Bromide	-	EPA Method 300.1	-	\$ 12.00
Chloride	-	EPA Method 300.1	-	\$ 12.00
Fluoride	-	EPA Method 300.1	-	\$ 12.00
Nitrate	-	EPA Method 300.1	-	\$ 12.00
Nitrite	-	EPA Method 300.1	-	\$ 12.00
Phosphate, ortho	-	EPA Method 300.1	-	\$ 12.00
Sulfate	-	EPA Method 300.1	-	\$ 12.00
Perchlorate	-	EPA Method 314.1	-	\$ 60.00

Physical Properties

pH	-	EPA Method 150.1	-	\$ 10.00
Specific Conductance	-	EPA Method 120.1	-	\$ 10.00
Turbidity	-	EPA Method 180.1	-	\$ 15.00
Biological Oxygen Demand (BOD)	-	EPA Method 405.1	-	\$ 40.00
Total Dissolved Solids (TDS)	-	EPA Method 160.1	-	\$ 20.00
Total Soluble Solids (TSS)	-	EPA Method 160.1	-	\$ 20.00